



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Four Post Lift

**Product Code :**  
LBNY-0002-640006



**Description :**

Four Post Lift

**Technical Specification :**

Widely used to lift different small and medium size vehicles with weight below 4.0t in garage and workshop. With the help of safety latches, these lifts offer mechanical safety protection.

Driven by hydraulic cylinder in stable lifting and lowering

Connected by four steel cables; forced synchronized movement of the lift in order to effectively prevent the sloping of the vehicles

Equipped with rolling jack, turntable and sides lip, suitable for wheel alignment

Lifting Capacity: 4 Ton

Lifting Time: 45 s

Lifting Height: 1915 mm

Overall Length: 5512 mm

Motor Power: 3Hp

Four Post Lift, Four Post Lift Equipments, Four Post Lift Tools, Four Post Lift Tool Kits, Four Post Lift Manufacturers, Four Post Lift Suppliers from India, China, Kenya

---

# Laboratory Instrument India

