



Email : sales@laboratoryinstrumentindia.com

Product Name :
Flame Spread Apparatus

Product Code :
LBNY-0005-1350002



Description :

Features:

Ignition source with fine-adjustment valve
Stainless steel combustion chamber
Timer
Specimen holder
Gas controls and safety flashback device.

Technical Specification :

Measuring principle: Horizontal flammability test
Bunsen burner tubes diameter: 9.5mm & 10mm supplied
Operating temperature: $23 \pm 5^{\circ}\text{C}$, non-condensing environment
Sample holder: U-shaped frame 50mm wide by 330mm long
Cabinet dimensions (interior):
Height: 356mm
Width: 381mm
Depth: 203mm.

Laboratory Instrument India

