



Email : sales@laboratoryinstrumentindia.com

Product Name :
Modular Air Flow Bench

Product Code :
LBNY-0005-1370002



Description :

Standard Features:

Supplied with a comprehensive user guide
Union directives

Available Experiment Modules:

Round Turbulent Jet
Boundary Layer
Flow Visualisation
Tapped Aerofoil
Bernoullis Equation
Drag Force
Flow Around a Bend
Jet Attachment

Recommended Ancillary:

Multi-Tube Manometer

Technical Specification :

Nett dimensions: 740 mm front to back x 950 mm wide x 1900 mm high and 125 kg

Approximate packed dimensions:

2.37 m³ and 150 kg

Operating Conditions

Storage temperature range: 25C to +55C (when packed for transport)

Operating temperature range: +5C to +40C

Operating relative humidity range: 80% at temperatures Operating environment: Laboratory.



Laboratory Instrument India