



Email : sales@laboratoryinstrumentindia.com

Product Name :
Universal Testing Machin Rotating Fatigue Machine

Product Code :
LBNY-0005-1380002



Description :

Salient Features:

Easy Changeability: Easy change from plain to threaded and screwed specimens.

Loading accuracy as high + 1%.

Autographic recorder: Simultaneous roll autographic recorder supplied as standard to enable study of the behaviors of materials.

Speeds: Straining at variable speeds to suit wide range of materials.

Facilities for tests: Motor-driven threaded columns for quick and convenient adjustment of lower cross head to facilitate rapid fixing of tear specimen.

Simple and Safe: Simple to operate. Robust construction. Chrome plated metal components. Wide range of standard and special attachments / accessories available.

Ideal Dia.: High reading accuracy due to large size ideal design of dial.

Large columns: Large effective clearance between columns enables testing of standard specimen as well as structures.

Technical Specification :

Universal Testing Machin Rotating Fatigue Machine

Laboratory Instrument India

