



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Structures Test Frame

**Product Code :**  
LBNY-0005-1390002



**Description :**

The frame has specially designed slots and self-positioning nuts that hold the Structures experiments and instruments. This fixing system is quick and easy to use. It allows students to change, position and secure each experiment.

A lightweight yet strong bench-mounting frame that holds interchangeable experiment modules and instrumentation from the All Structures Tester range.

Adjustable feet support the frame to allow students to level the apparatus before use.

Supplied in kit form with instructions for use and a text book.

**Technical Specification :**

**Features:**

Easy-to-use fixings and slots so students can quickly set up, remove or change experiments

Also ideal for holding experiments during storage

A strong frame that holds the experiments of Alls Structures Test range

Strong, bench-mounting frame

Holds the interchangeable experiment modules and instruments of Alls Structures Test range.



**Laboratory Instrument India**

