



Email : sales@laboratoryinstrumentindia.com

Product Name :
Pendulum Impact Tester

Product Code :
LBNY-0002-640008



Description :

Pendulum Impact Tester

Technical Specification :

Semi-automatic working and high energy capacity
Mechanical safety with automatic insertion to the arm hooking.
Simple and fast utilization, ideal for routine tests.
Supplied complete with protection cage to CE Safety Directive.
Impact energy: 300J with 0.25J resolution.
Power supply: 230V 1hp,50 Hz 180W
Dimensions: 800 x578x 1400mm
Weight: 450kg.

Pendulum Impact Tester, Pendulum Impact Tester Equipments, Pendulum Impact Tester Tools, Pendulum Impact Tester Tool Kits, Pendulum Impact Tester Manufacturers, Pendulum Impact Tester Suppliers from India, China, Kenya



Laboratory Instrument India

