



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Pendulum Impact Tester

**Product Code :**  
LBNY-0002-640008



**Description :**

Pendulum Impact Tester

**Technical Specification :**

Semi-automatic working and high energy capacity  
Mechanical safety with automatic insertion to the arm hooking.  
Simple and fast utilization, ideal for routine tests.  
Supplied complete with protection cage to CE Safety Directive.  
Impact energy: 300J with 0.25J resolution.  
Power supply: 230V 1hp,50 Hz 180W  
Dimensions: 800 x578x 1400mm  
Weight: 450kg.

Pendulum Impact Tester, Pendulum Impact Tester Equipments, Pendulum Impact Tester Tools, Pendulum Impact Tester Tool Kits, Pendulum Impact Tester Manufacturers, Pendulum Impact Tester Suppliers from India, China, Kenya



**Laboratory Instrument India**

