



Email : sales@laboratoryinstrumentindia.com

Product Name :
Special Weight Length Minus Tolerance Measuring Test Equipment

Product Code :
LBNY-0005-1450001



Description :

Length measurement accuracy 0.2mm Special Weight Length Minus Tolerance Measuring Test Equipment With High Precision.

Features:

Specially developed measurement control system, featuring a high speed processing unit and a performance of high stability and precision.
Safety barriers are installed at two sides of the working table to ensure the safety of operators.
Integrated with non-contact laser measurement, improving the precision and reliability of measurement data.
Completely automatic measurement process. Place the component on the V-shaped bearing, press the measure button, and the measurement instrument will complete the measurement quickly. The measurement data are directly retrieved, displayed and saved by computer software.
Integrated with safety barriers, travel limit protection and full scale protection.
Intuitive display of length and weight.
Support for saving, instant retrieving and inspecting measurement data.

Technical Specification :

Special Weight Length Minus Tolerance Measuring Test Equipment

Laboratory Instrument India

