



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Special Weight Length Minus Tolerance Measuring Test Equipment

**Product Code :**  
LBNY-0005-1450001



**Description :**

Length measurement accuracy 0.2mm Special Weight Length Minus Tolerance Measuring Test Equipment With High Precision.

**Features:**

- Specially developed measurement control system, featuring a high speed processing unit and a performance of high stability and precision.
- Safety barriers are installed at two sides of the working table to ensure the safety of operators.
- Integrated with non-contact laser measurement, improving the precision and reliability of measurement data.
- Completely automatic measurement process. Place the component on the V-shaped bearing, press the measure button, and the measurement instrument will complete the measurement quickly. The measurement data are directly retrieved, displayed and saved by computer software.
- Integrated with safety barriers, travel limit protection and full scale protection.
- Intuitive display of length and weight.
- Support for saving, instant retrieving and inspecting measurement data.

**Technical Specification :**

Special Weight Length Minus Tolerance Measuring Test Equipment

---

# Laboratory Instrument India

