



Email : sales@laboratoryinstrumentindia.com

Product Name :
Hydraulic Pressure To Clamp For Sheet Metal Ductility
Testing Machine

Product Code :
LBNY-0005-1460002



www.laboratoryinstrumentindia.com

Description :

Cup drawing Display error 0.01mm
Max force 300kN
Clamp force 300kN
Clamping force relative error 2%
Test force range 2%100%
Test force display error 1%
Displacement resolution 0.002mm
Stamping stroke 120mm
Clamp stroke 40mm
Stamping speed 130mm/min
Plate maximum width 200mm
power 3.7kW/AC380V/50Hz
frame 930mmx1260mmx1400mm
weight About 1800kg

Technical Specification :

Hydraulic Pressure To Clamp For Sheet Metal Ductility Testing Machine

Laboratory Instrument India

