



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Compression Testing Machine Frame Type With Hydraulic

**Product Code :**  
LBNY-0005-1480001



**Description :**

Domestically pioneered open-testing-method settings. Engineers may set up testing projects then operators may employ the projects for testing, greatly reduced the difficulty of the operations and the requirement in the skills and trainings of the operator.

Domestically pioneered highly efficient three-level operation model of testing data processing. Supports open testing data processing, appending custom testing data items and re-inspection of processed testing data. Dynamic display of load, displacement, time and real-time testing curve on the screen.

Support for compression tests on various metal or non-metal materials to inspect their mechanics performances.

Domestically pioneered process-oriented design of testing control. Provides users with intuitive testing control procedures and a wide range of instruction types. Also supports dynamic instruction parameter assignment.

Closed loop control of channels such as load, displacement, etc. Smooth switching in between.

Constant loading speed test supported.

Capable of holding on maximum load or displacement.

Characteristic testing data can be marked on testing graphs.

Capable of statistical operations on testing data .

Automatic requesting for testing result data, also capable of launching an advanced analysis to inspect the data.

**Technical Specification :**

Compression Testing Machine Frame Type With Hydraulic

---

# Laboratory Instrument India

