



Email : sales@laboratoryinstrumentindia.com

Product Name :
Hydraulic Servo Compression Testing Machine for Testing

Product Code :
LBNY-0005-1480002



Description :

Features:

Dynamic display of load, displacement, time and real-time testing curve on the screen.
Constant loading speed test supported.
Capable of holding on maximum load or displacement.
Support for multi-language user-interface extendibility.
Closed loop control of channels such as load, displacement, etc. Smooth switching in between.

Technical Specification :

Hydraulic Servo Compression Testing Machine for Testing



Laboratory Instrument India