



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Hydraulic Servo Compression Testing Machine for Testing

**Product Code :**  
LBNY-0005-1480002



**Description :**

**Features:**

Dynamic display of load, displacement, time and real-time testing curve on the screen.  
Constant loading speed test supported.  
Capable of holding on maximum load or displacement.  
Support for multi-language user-interface extensibility.  
Closed loop control of channels such as load, displacement, etc. Smooth switching in between.

**Technical Specification :**

Hydraulic Servo Compression Testing Machine for Testing



**Laboratory Instrument India**