



Email : sales@laboratoryinstrumentindia.com

Product Name :
Combustion Laboratory Unit

Product Code :
LBNY-0005-1500001



www.laboratoryinstrumentindia.com

Description :

Internal electric and mechanical safety devices allow for supervised operation by students. Instrumentation to measure the water, air and exhaust temperatures, and flow rates of the cooling water, air and fuel. A hand held gas analyser supplied measures flue gas O₂ content and also displays combustion efficiency, excess air, and CO₂ content.

A frame mounted burner training unit comprising a 450 dia x 1000m cylindrical water-cooled stainless steel chamber in which either an oil or gas burner can be fired, and the flame observed through four 100mm dia windows.

Oil or Gas Burners up to 150kW can be used in the chamber.

Technical Specification :

Supplied with a detailed experimental operating and maintenance manual giving example experimental results and sample calculations.

Smoke Number can be determined for the oil burner with the Smoke Tester supplied.

Further flue gas analysis including CO, NO and SO₂ can be performed on either burner by using the additional optional high specification Gas Analyser.



Laboratory Instrument India

