



Email : sales@laboratoryinstrumentindia.com

Product Name :
Flame Propagation And Stability Unit

Product Code :
LBNY-0005-1500002



www.laboratoryinstrumentindia.com

Description :

Panel: vacuum formed solid colour panel in which all controls, fan and instruments are mounted

Burner Tubes: 4 inter-changeable steel burner tubes in a range of diameters.

Stabiliser Cups: two brass conical cups to fit two of the burner tubes.

Airflow: Glass variable area flowmeter

Gas Flow: Glass variable area flowmeter to cover wide range of load speed combustible gases(e.g. Propane, Butane, Methane LPG).

Flame Speed System: Long, large diameter clear tube, with mains driven high-voltage spark igniter and two adapters each with flame traps.

Technical Specification :

Smithells Tube: heat resistant glass tube and brass adapter.

Safety: Gas control solenoid valve with foot switch. Combined overload cut out and on /off switch, earth leakage circuit breaker.

Igniter: Handheld manual.



Laboratory Instrument India

