



Email : sales@laboratoryinstrumentindia.com

Product Name :
Hydraulic Fatigue Testing Machine With Dynamic Display

Product Code :
LBNY-0005-1510002



Description :

The precise actuator collaborates with the intelligent Test Live computer testing software to ensure the completeness and precision of the test results. This product is integrated with a high-stiffness frame. The actuator cylinder comes out in 100mm, 150mm, 400mm and many other models. The crossbeam is hydraulically driven to adjust the test space. The rated load of the main frame ranges from 5 to 25kN or 50 to 1000kN.

Technical Specification :

Hydraulic Fatigue Testing Machine With Dynamic Display



Laboratory Instrument India