



Email : sales@laboratoryinstrumentindia.com

Product Name :
Multifunction Fastening Analysis System Testing Machine
Unit

Product Code :
LBNY-0005-1520001



Description :

Max capacity 1000 Nm
Max thread of friction torque 800 Nm
Torque measurement range 20 Nm ~1000 Nm
Torque accuracy $\pm 1\%$
Torque resolution ± 500000 yard, whole resolution unchanged
Max clamp force 300 kN
Clamp force accuracy $\pm 0.5\%$
Clamp force measurement range 3 kN ~300 kN
Clamp force resolution ± 500000 yard, whole resolution unchanged
Torsional angle accuracy $\pm 1\%$
Torsional angle measurement resolution 0.025°
Max torsion speed 15 r/min
Torsion speed accuracy $\pm 1\%$
Torsion direction Positive and negative two directions
Max test length of fasteners 300mm
Power supply 2.5Kw/ AC220v/50Hz
Frame dimensions 2000mm \times 600mm \times 1200mm
Weight 860kg

Technical Specification :

Multifunction Fastening Analysis System Testing Machine Unit



Laboratory Instrument India