



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Computer Controlled Electronic Torsion Testing Machine

**Product Code :**  
LBNY-0005-1520002



**Description :**

Automatically sustains when it reaches the specified value.  
Capable of re-analyzing test data and curves, and partial zooming of testing curves at will.  
Support Excel report templates and customized editing of test reports.  
Support for data searching. Returns the corresponding data at y-axis when entered the coordinate at x-axis.  
Support saving, instant retrieving and inspection of testing results.  
Dynamic display of torque, rotation and real-time testing curve on the screen.  
Integrated with travel limit protection and full-scale overload protection.  
Capable of sustaining at the maximum torque, rotation .  
Options for displaying curves including torque-rotation.  
Support for setting a target value on torque, rotation and grip load at will.  
Control over loading speed of torque and rotation, no stepping for loading speed.

**Technical Specification :**

Computer Controlled Electronic Torsion Testing Machine



**Laboratory Instrument India**

