



Email : sales@laboratoryinstrumentindia.com

Product Name :
Space Rebar Cold Bending Test Machine

Product Code :
LBNY-0005-1530001



Description :

Safe: specimen bending directly into sealed test table, safe to operate.
Dust-proof: unique dust resistant unit protects ball screw from oxide coating.
Two-way hydraulic cylinder enables fast loading, unloading and withdrawing.
Effortless: requiring only one operator, who clicks the buttons to perform tests.
Equipped with waste collecting cart, reduces working effort.
Small footprint saves space for testing.
One-step forming: this product bends specimen in one step during tests. Then applies side loads to form interference bend on the specimen to eliminate remaining deformation, accomplishing a bending angle of 180.

Technical Specification :

Space Rebar Cold Bending Test Machine



Laboratory Instrument India