



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Bow Type Drop Optical Cable Tensile Testing Machine

**Product Code :**  
LBNY-0005-1570001



**Description :**

This product features a brand new structural design, optimizes and upgrades on the base of last generation. Remote communication ports are installed to make it easier for remote controlling operation on a computer. Impact height may be configured at will, fulfilling all requirements for optical cable impact tests under 1000 mm.

**Technical Specification :**

Impact count: 19999 times, configured at will.  
Impact speed: 1030 times/min, configured at will.  
Hammer weight: 16kg in all. Weights: 0.5kg, 2, 1kg, 3, 3kg, 2, 6kg, 1.  
Tup insert with bracket weight: 0.55kg.



**Laboratory Instrument India**