



Email : sales@laboratoryinstrumentindia.com

Product Name :
Bow Type Drop Optical Cable Tensile Testing Machine

Product Code :
LBNY-0005-1570001



Description :

This product features a brand new structural design, optimizes and upgrades on the base of last generation. Remote communication ports are installed to make it easier for remote controlling operation on a computer. Impact height may be configured at will, fulfilling all requirements for optical cable impact tests under 1000 mm.

Technical Specification :

Impact count: 19999 times, configured at will.
Impact speed: 1030 times/min, configured at will.
Hammer weight: 16kg in all. Weights: 0.5kg, 2, 1kg, 3, 3kg, 2, 6kg, 1.
Tup insert with bracket weight: 0.55kg.



Laboratory Instrument India