



Email : sales@laboratoryinstrumentindia.com

Product Name :
Water Pumps Assembly Testing Bench

Product Code :
LBNY-0005-1600002



Description :

The Water Pumps Assembly Testing Bench teaching unit is a modular group designed to perform complete functional tests on a series of 4 water pumps. Each pump is taken into consideration individually and with it the following characteristic data can be obtained:

Suction pressure
Power consumption
Delivery pressure
Speed
Flow rate

Technical Specification :

Weight and Dimensions:
Control board: 350x500x400 h mm
Weight: 300 kg.



Laboratory Instrument India

