



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Water Pumps Assembly Testing Bench

**Product Code :**  
LBNY-0005-1600002



**Description :**

The Water Pumps Assembly Testing Bench teaching unit is a modular group designed to perform complete functional tests on a series of 4 water pumps. Each pump is taken into consideration individually and with it the following characteristic data can be obtained:

Suction pressure  
Power consumption  
Delivery pressure  
Speed  
Flow rate

**Technical Specification :**

Weight and Dimensions:  
Control board: 350x500x400 h mm  
Weight: 300 kg.



**Laboratory Instrument India**

