



Email : sales@laboratoryinstrumentindia.com

Product Name :
Technique Training Extend PLC Model

Product Code :
LBNY-0002-650004



Description :

Technique Training Extend PLC Model

Technical Specification :

This training device can conduct PLC typical application project training.

Dimension: About 1500 x 700 x 1700mm

Working power: AC 220V

AC/DC Power Supply: 0-24v(DC)/220V, 1A AC

Technique Training Extend PLC Model, Technique Training Extend PLC Model Equipments, Technique Training Extend PLC Model Tools, Technique Training Extend PLC Model Tool Kits, Technique Training Extend PLC Model Manufacturers, Technique Training Extend PLC Model Suppliers from India, China, Kenya



Laboratory Instrument India