



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Auto-Trainer Engine Test Stands

**Product Code :**  
LBNY-0005-1690002



[www.laboratoryinstrumentindia.com](http://www.laboratoryinstrumentindia.com)

**Description :**

The test stands are supplied complete with cooling, electrical, exhaust and fuel systems with provision for the optional factory fitted Instrumentation. Auto-trainer engine test stands are of all welded construction, made from square section steel tube, and are fully mobile being mounted on four castors, two of which are lockable. The engines are flexibly mounted to reduce vibration and guarded as necessary.

**Applications:**

Engine Operation, Maintenance and Testing  
Fault Diagnostics  
Idle Speed Adjustment

**Features:**

Naturally Aspirated or Turbo-charged Engines  
Petrol, LPG or Diesel Engine Variants  
Single, in line and V cylinder configuration.

**Technical Specification :**

Standard Instruments And Controls of Auto-Trainer Engine Test Stands:  
The unit incorporates the necessary instruments, which may include some of the following :  
Water temperature gauge,  
Pre-start indicator - Diesel,  
Manual throttle or rack control,  
Tachometer and ammeter,

---

Ignition-starter key switch,  
Alternator charge warning light,  
Oil pressure gauge.



**Laboratory Instrument India**