



Email : sales@laboratoryinstrumentindia.com

Product Name :
Industrial Multimeter

Product Code :
LBNY-0002-660001



Description :

Industrial Multimeter, Industrial Multimeter Equipments, Industrial Multimeter Tools, Industrial Multimeter Tool Kits, Industrial Multimeter Manufacturers, Industrial Multimeter Suppliers from India, China, Kenya

Lab Tender Suppliers Africa and Lab Equipment Africa Kenya, Zambia, Chad, Somalia, Senegal, Zimbabwe, Guinea, Rwanda, Benin, Burundi, Tunisia, Uganda, Sudan, Algeria, Morocco, Angola, Ghana, Mozambique, Madagascar, South Africa, Ivory Coast, Cameroon, Niger, Burkina Faso, Mayotte, Sao Tome, Seychelles, South Sudan, Mali, Malawi, Nigeria, Ethiopia, Egypt, DR Congo, Mauritania, Eritrea, Gambia, Botswana, Namibia, Gabon, Lesotho, Guinea Bissau, Equatorial Guinea, Mauritius, Eswatini, Djibouti, Reunion, Comoros, Cape Verde, Western Sahara, Tanzania, Togo, Sierra Leone, Libya, Congo, Central African Republic, Liberia

Technical Specification :

4000 Counts HD Display with backlight, True RMS Auto / Manual ranging and Shockproof construction which ensures ultimate protection

Duty-Cycle, Frequency in ACV range, Relative measurements, Transistor (hFE) test, Capacitance 6 nF to 100 mF Temperature - 20 °C ~ 1000 °C Resistance 400 Ω to 40 M Ω Frequency 10 Hz to 20 MHz

High stability high precision instrument with full function protection,

Low battery indicator, auto power off, data hold

Power Source: Battery Powered

AC/DC: 1000 V

Style: Adjustable, Unique, Magnetic

Colour: Industrial True RMS Multimeter.



Laboratory Instrument India