



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Laminar Flow Meter Range Instrument

**Product Code :**  
LBNY-0005-1700001



[www.laboratoryinstrumentindia.com](http://www.laboratoryinstrumentindia.com)

**Description :**

**Experimental Ability:**

Single or multi-cylinder engines  
Naturally aspirated or supercharged engines  
Measure the air flow of the engine, hose, or leads  
Air Compressors  
Research and development.

**Technical Specification :**

**Functionalities:**

Can be used with pulsating flows  
Not affected by elbows or bends in conduits  
Verified Linear Characteristics  
Accurate measurement of flow  
Ranges available from 25 to 650 liters per second  
Improved pneumatic damping at connections Pressure  
"Square Law" Errors eliminated by adjustment  
Removable element for cleaning.

---

# Laboratory Instrument India

