



Email : sales@laboratoryinstrumentindia.com

Product Name :
Exhaust Gas Analyser

Product Code :
LBNY-0005-1710002



Description :

Warm up time 30 seconds
Bluetooth, wireless connection for PC/Notebook or handheld as an option
Internal OBD-scantool as an option
Small, robust and compact gas analyser designed for workshop environments or mobile use
Available with an easy to use handset or PC based program Easy to use, low user maintenance

Technical Specification :

Measuring ranges:
CO₂ 0 - 20 % Vol.
CO 0 - 10 % Vol.
O₂ 0 - 22 % Vol.
HC 0 - 2500 ppm Vol.
Working temperature: 5 - 45°C
Warm up time: 30s
Lambda calculation: 0 - 9,999
Temperature measuring range: 0 - 130 °C
RPM measuring range: 400 - 9999 min⁻¹.

Laboratory Instrument India

