



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Insulation Multimeter

**Product Code :**  
LBNY-0002-660002



#### **Description :**

Insulation Multimeter

#### **Technical Specification :**

A versatile, compact, handheld insulation tester for advanced industrial and electrical insulation testing. Its multiple test voltages make it the ideal tool for many troubleshooting, commissioning, and preventative maintenance applications.

Take advantage of the full feature set to easily and quickly perform advanced insulation resistance testing. Its handy size and light weight make it easy to pack and use.

Additional features, like the remote probe, save both time and money especially when performing repetitive tests.

And its reasonable price and rugged reliability make it an excellent value.

Power Source: Corded Electric

Insulation range : 0.01 M $\Omega$  to 2 G $\Omega$

Insulation test voltage: 50, 100, 250, 500, 1000 V.

Insulation Multimeter, Insulation Multimeter Equipments, Insulation Multimeter Tools, Insulation Multimeter Tool Kits, Insulation Multimeter Manufacturers, Insulation Multimeter Suppliers from India, China, Kenya

---

# Laboratory Instrument India

