



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Beaker Low Form

**Product Code :**  
LBNY-0005-1750001



**Description :**

Capacity (ml) D. (mm) H. (mm)

5 22 30

10 26 35

25 34 50

50 42 60

100 50 70

100 60 80

250 70 95

400 80 110

600 90 125

1000 105 145

2000 132 185

**Technical Specification :**

Beaker Low Form



**Laboratory Instrument India**

