



Email : sales@laboratoryinstrumentindia.com

Product Name :
Linear Motion Test Apparatus

Product Code :
LBNY-0005-1770001



Description :

Bench top for investigations on motion, plastic and elastic impact, momentum, spring-mass oscillation and on friction System to consist of:

Inclinable track with stoppers at both sides

Friction block

Plunger carriage

Device for actuation of carriages and friction block

3 tension springs

Rolling carriage

3 sets of weights

Deflection pulley

Cord

Moveable stopper for track

3 buffers with rubber buffer heads for covering of permanent magnets

Additional weights for carriages including measuring tape: length: 2000mm, resolution: 1mm.

Technical Specification :

Linear Motion Test Apparatus

Laboratory Instrument India

