



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Dynamometers For Educational Towing Carriages

**Product Code :**  
LBNY-0005-1780002



[www.laboratoryinstrumentindia.com](http://www.laboratoryinstrumentindia.com)

**Description :**

A variety of dynamometers can be supplied for use with the educational carriages, these are interfaced with the support frame for connection to ship models or fitting of propeller models. These include:

- Educational Multi Component Ship / Yacht Model Balance
- Educational Ship Stability Dynamometer
- Educational Self Propulsion
- Traditional Self-Propulsion Dynamometer for educational scaled models
- Educational Resistance Dynamometer
- Educational Open Water Propeller / Turbine Dynamometer.

**Technical Specification :**

Dynamometers For Educational Towing Carriages



**Laboratory Instrument India**