



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Insulated Offset End Wrench Set (15 PCS)

**Product Code :**  
LBNY-0002-660005



**Description :**

Insulated Offset End Wrench Set (15 PCS)

**Technical Specification :**

Dual color nonslip handle for max. comfort

Made of Chrome Vanadium Steel.

With 600D tool roll bag.

Ideal for power plant, HVAC, automobile, military, aircraft field applications.

(8mm, 9mm, 10mm, 11mm, 12mm, 13mm, 14mm, 15mm, 16mm, 17mm, 18mm, 19mm, 21mm, 22mm, 24mm).

Insulated Offset End Wrench Set (15 PCS), Insulated Offset End Wrench Set (15 PCS) Equipments, Insulated Offset End Wrench Set (15 PCS) Tools, Insulated Offset End Wrench Set (15 PCS) Tool Kits, Insulated Offset End Wrench Set (15 PCS) Manufacturers, Insulated Offset End Wrench Set (15 PCS) Suppliers from India, China, Kenya



**Laboratory Instrument India**

