



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Oven 440 Liters Natural Convection

**Product Code :**  
LBNY-0005-1950001



**Description :**

The oven is supplied complete with two grid shelves easily removable and that can be positioned at various heights pilot light, exhaust holes for fast cooling. Natural convection digital thermoregulator indicator. Designed for drying, baking, conditioning and moisture determination. Sturdy manufacture double walled with 60 mm thick glass fiber for thermal insulation. Exterior front part is stainless steel made; while interior chamber, grid shelves and external walls are made from zinc coated steel. The oven is equipped of a dual safety thermostat with higher thermic threshold to prevent accidental over-temperatures, and to ensure safe working conditions. Temperature from ambient to 200°C. is controlled by a digital thermoregulator-indicator.

**Technical Specification :**

Power supply: 230V 50-60Hz 1ph  
Wattage: 3600W  
Doors: n.2  
Inside dimensions: 900X700X700 mm  
Outside dimensions: 1140X760X910 mma  
Weight: 85 kg.



**Laboratory Instrument India**

