



Email : sales@laboratoryinstrumentindia.com

Product Name :
Oven 220 Liters Natural Convection

Product Code :
LBNY-0005-1950002



Description :

Designed for drying, baking, conditioning and moisture determination. Sturdy manufacture, double walled with 60 mm thick glass fiber for thermal insulation. Exterior front part is stainless steel made; while interior chamber, grid shelves and external walls are made from zinc coated steel. Natural convection, digital thermoregulator-indicator. Temperature from ambient to 200°C. is controlled by a digital thermoregulator-indicator. The oven is supplied complete with two grid shelves easily removable and that can be positioned at various heights, pilot light, and exhaust holes for fast cooling. The oven is equipped of a dual safety thermostat with higher thermic threshold to prevent accidental over-temperatures, and to ensure safe working conditions.

Technical Specification :

Doors: n.1
Wattage: 2000W
Power supply: 230V 50-60Hz 1ph
Inside dimensions: 600X610X600 mm
Outside dimensions: 840X725X805 mm
Weight: 60 kg.



Laboratory Instrument India

