




Email : sales@laboratoryinstrumentindia.com

Product Name :
Function Generator

Product Code :
LBNY-0002-670001



Description :

Function Generator, Function Generator Equipments, Function Generator Tools, Function Generator Tool Kits, Function Generator Manufacturers, Function Generator Suppliers from India, China, Kenya

Lab Tender Suppliers Africa and Lab Equipment Africa Kenya, Zambia, Chad, Somalia, Senegal, Zimbabwe, Guinea, Rwanda, Benin, Burundi, Tunisia, Uganda, Sudan, Algeria, Morocco, Angola, Ghana, Mozambique, Madagascar, South Africa, Ivory Coast, Cameroon, Niger, Burkina Faso, Mayotte, Sao Tome, Seychelles, South Sudan, Mali, Malawi, Nigeria, Ethiopia, Egypt, DR Congo, Mauritania, Eritrea, Gambia, Botswana, Namibia, Gabon, Lesotho, Guinea Bissau, Equatorial Guinea, Mauritius, Eswatini, Djibouti, Reunion, Comoros, Cape Verde, Western Sahara, Tanzania, Togo, Sierra Leone, Libya, Congo, Central African Republic, Liberia

Technical Specification :

Frequency response analyzer & function generator dual function
Wide sweep frequency up to 60MHz
Dual channels function/arbitrary waveform generation
Channel independence, coupling, track working mode
Complete dual channels function/arbitrary waveform generator
Channel independence, coupling, track working mode
200 MSa/s sampling rate and 14-bit vertical resolution per channel
Output of 6 standard waveforms, built-in 50 kinds of arbitrary waveform
Multi modulation function: AM, DSSC - AM, FM, PM, PWM, FSK, ASK, BPSK and logarithm/linear sweep.



Laboratory Instrument India