



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Tricycle Landing Gear Trainer

**Product Code :**  
LBNY-0005-2000002



**Description :**

Tricycle Landing Gear Trainer can be used by students for hands-on training related to light aircraft landing gears.

**The students task may include:**

Inspection and maintenance of nose wheel steering mechanism

The Tricycle Landing Gear Trainer is fitted with actual aircraft parts that includes tires, tubes, wheels, dual Cleve and single-disc brakes, nose-wheel steering mechanism, brake pedals, dual brake-fluid reservoirs, master cylinders and hydraulic plumbing for the brakes.

Shock strut inspection and maintenance

One of the three shock struts has a cutaway to expose the internal mechanism, while the other two are complete, and functional shock struts.

Shock strut is assembly and reassembly, and Hydraulic brakes inspection and maintenance.

Dimensions: 48L x 48W x 64H inches.

**Technical Specification :**

Tricycle Landing Gear Trainer

---

# Laboratory Instrument India

