



Email : sales@laboratoryinstrumentindia.com

Product Name :
Oscilloscope (100 MHz)

Product Code :
LBNY-0002-670002



Description :

Oscilloscope (100 MHz)

Technical Specification :

Automatic measurement of 34 waveform parameters
Automatic cursor tracking measurement function
2 channels, low noise floor, wide vertical range: 1mV/div~20V/div
Built in FFT and digital filter
System software upgrade via USB drive
7 inches TFT LCD,high resolution display 800×480
Delicate window extension function and precise analysis on waveform details and overview
AUTO's strategy presetable
Automatic measurement of 34 waveform
Automatic cursor tracking measurement
Low noise floor, wide vertical range:
7 inches TFT LCD,high resolution display.

Oscilloscope (100 MHz), Oscilloscope (100 MHz) Equipments, Oscilloscope (100 MHz) Tools, Oscilloscope (100 MHz) Tool Kits, Oscilloscope (100 MHz) Manufacturers, Oscilloscope (100 MHz) Suppliers from India, China, Kenya



Laboratory Instrument India