




Email : sales@laboratoryinstrumentindia.com

Product Name :
Oscilloscope (100 MHz)

Product Code :
LBNY-0002-670002



Description :

Oscilloscope (100 MHz)

Technical Specification :

Automatic measurement of 34 waveform parameters

Automatic cursor tracking measurement function

2 channels, low noise floor, wide vertical range: 1mV/div~20V/div

Built in FFT and digital filter

System software upgrade via USB drive

7 inches TFT LCD,high resolution display 800x480

Delicate window extension function and precise analysis on waveform details and overview

AUTO's strategy presetable

Automatic measurement of 34 waveform

Automatic cursor tracking measurement

Low noise floor, wide vertical range:

7 inches TFT LCD,high resolution display.

Oscilloscope (100 MHz), Oscilloscope (100 MHz) Equipments, Oscilloscope (100 MHz) Tools, Oscilloscope (100 MHz) Tool Kits, Oscilloscope (100 MHz) Manufacturers, Oscilloscope (100 MHz) Suppliers from India, China, Kenya



Laboratory Instrument India