



Email : sales@laboratoryinstrumentindia.com

Product Name :
Avionics System Trainer

Product Code :
LBNY-0005-2010002



Description :

The trainer is designed to give the students a good understanding aircrafts avionics equipment. Students also develop logical and methodical approach to trouble shoot avionics faults inserted by the instructor. Avionics System Trainer brings the real life experience to the classroom. The trainer provides practical, hands-on training using the real avionics equipment configured and wired just like in an airworthy aircraft.

Standard Equipment:

Air craft avionics trainer provides instruction on the following systems

- Global Positioning System
- Altitude Encoder
- Automatic Direction Finder
- Instrument Landing System
- Dual NAV. /Comm. Radio
- Transponder
- DME
- Marker Beacon
- Microphone
- Two NAV Indicators
- ADF Indicator
- VOR
- Intercom System
- Headsets
- Audio Panel
- Fault Insertion Box

LCD NAV Display with GPS built in
Altimeter
Antennas.

Technical Specification :

Avionics System Trainer



Laboratory Instrument India