



Email : sales@laboratoryinstrumentindia.com

Product Name :
Consolidation Apparatus

Product Code :
LBNY-0005-2020002



Description :

The load bridge group is supported in high accuracy self-aligning seat balls. Rigidly manufactured from aluminous alloy casting to provide a high degree of accuracy with any frame distortions under load.

Technical Specification :

It consists of 3 consolidation test assemblies mounted on a single frame. Complete with 3 sets of Consolidation cell assembly for 60mmdia samples and three sets of weights each giving 10 kg/sq.cm. The beam provides three load ingratio: 9:1 10:1 11:1 and the beam assembly is fitted with an adjustable counterbalance weight. On each sample And complete with rod holding the weights and coupling block holding the dial gauge or transducer.



Laboratory Instrument India