



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
FMS Mechatronics Training System

**Product Code :**  
LBNY-0006-2060001



**Description :**

FMS Mechatronics Training System, FMS Mechatronics Training System Equipments, FMS Mechatronics Training System Tools, FMS Mechatronics Training System Tool Kits, FMS Mechatronics Training System Manufacturers, FMS Mechatronics Training System Suppliers from India, China, Kenya

Lab Tender Suppliers Africa and Lab Equipment Africa Kenya, Tunisia, Uganda, Sudan, Algeria, Morocco, Angola, Ghana, Gabon, Lesotho, Zambia, Chad, Somalia, Senegal, Zimbabwe, Guinea, Rwanda, Benin, Burundi, Guinea Bissau, Equatorial Guinea, Mauritius, Eswatini, Djibouti, Reunion, Comoros, Cape Verde, Western Sahara, Tanzania, Togo, Sierra Leone, Libya, Congo, Central African Republic, Liberia, Mozambique, Madagascar, South Africa, Ivory Coast, Cameroon, Niger, Burkina Faso, Mayotte, Sao Tome, Seychelles, South Sudan, Mali, Malawi, Nigeria, Ethiopia, Egypt, DR Congo, Mauritania, Eritrea, Gambia, Botswana, Namibia

**Technical Specification :**

Main performance characteristics Use 220 volts, fully enclosed plus transparent plexiglass structure, high-quality cast iron material casting, high-precision grinding ball screw;  
Equipped with Yornew980TB industrial panel CNC system;  
Execute the international standard G code programming, support M code and S code, compatible with FANUC, Mitsubishi G code and various CAD / CAM software (MasterCAM, UG, CAXA, etc. software programming, etc.);

---

# Laboratory Instrument India

