



Email : sales@laboratoryinstrumentindia.com

Product Name :
Portable Cutting Machine

Product Code :
LBNY-0002-670003



Description :

Portable Cutting Machine

Technical Specification :

Less fatigue: Compression spring design saves energy during long operation hours.

Easy to carry: Reduced tool weight by 1-2kg with optimized base design.

Precise cutting: Optimized guard designs reduces tools vibration and optimized base design leads to improved cutting precision.

Portable: Compact size, portable and easy to carry.

Power: 1500W.

Voltage: 220V.

Speed: 11000 rpm.

Saw Type: Mini Cutting Machine.

Portable Cutting Machine, Portable Cutting Machine Equipments, Portable Cutting Machine Tools, Portable Cutting Machine Tool Kits, Portable Cutting Machine Manufacturers, Portable Cutting Machine Suppliers from India, China, Kenya

Laboratory Instrument India

