



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Manipulator Training Equipment

**Product Code :**  
LBNY-0006-2060002



**Description :**

Manipulator Training Equipment, Didactic Equipment, Educational

**Technical Specification :**

This trainer takes various kinds of industrial sensor, pneumatic controller, electrical switch, PLC, HMI etc., it is a comprehensive mechatronics platform which consist of warehouse unit, handling mechanical arm(manipulator) unit, movement control PLC module, touch screen unit etc. Through related experiment, students will be familiar with pneumatic execute mechanism running feature, Siemens 1200PLC based movement control system, touch screen and PLC communication method etc. , to cultivate students corresponding knowledge and technical skills, it is suitable for vocational schools, colleges, technical institutions teaching and training for hand on practice ability.



**Laboratory Instrument India**