



Email : sales@laboratoryinstrumentindia.com

Product Name :
General Air Conditioning Trainer

Product Code :
LBNY-0006-2070002



Description :

General Air Conditioning Trainer, HVAC Training Equipment

Technical Specification :

The training device is designed for the in-depth study of thermodynamics, make the air go through the different stages of air conditioning refrigeration unit, and then enter a room which temperature and humidity can be real-time monitoring, measuring the difference of air temperature and humidity to analysis of air cooling, heating, humidification and dehumidification, to make the students to understand the feature of open compressor, study the specific operation and the working principle of an open compressor. In addition, some indoor simulator and outdoor temperature contains the logic control system in all of the environment, can occur in air handling unit.



Laboratory Instrument India