



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Didactic Equipment Hydrostatic Trainer

**Product Code :**  
LBNY-0006-2080002



**Description :**

Didactic Equipment Hydrostatic Trainer

**Technical Specification :**

Power consumption: approx. 250W  
max. flow rate: 9m<sup>3</sup>/h  
max. head: 7,6m  
Compressor  
power: approx. 65W  
delivery side: 2000mbar  
intake side: 240mbar  
tanks  
Supply tank for water: approx. 55L.  
areometers with different measuring ranges



**Laboratory Instrument India**

