



Email : sales@laboratoryinstrumentindia.com

Product Name :
Didactic Equipment Hydrostatic Trainer

Product Code :
LBNY-0006-2080002



Description :

Didactic Equipment Hydrostatic Trainer

Technical Specification :

Power consumption: approx. 250W
max. flow rate: 9m³/h
max. head: 7,6m
Compressor
power: approx. 65W
delivery side: 2000mbar
intake side: 240mbar
tanks
Supply tank for water: approx. 55L.
areometers with different measuring ranges



Laboratory Instrument India

