




Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Speed and Position Transducer and Control Trainer

**Product Code :**  
LBNY-0006-210001



#### **Description :**

Speed and Position Transducer and Control Trainer, Speed and Position Transducer and Control Trainer Equipments, Speed and Position Transducer and Control Trainer Tools, Speed and Position Transducer and Control Trainer Tool Kits, Speed and Position Transducer and Control Trainer Manufacturers, Speed and Position Transducer and Control Trainer Suppliers from India, China, Kenya

Lab Tender Suppliers Africa and Lab Equipment Africa Kenya, Tunisia, Uganda, Sudan, Algeria, Morocco, Angola, Ghana, Gabon, Lesotho, Zambia, Chad, Somalia, Senegal, Zimbabwe, Guinea, Rwanda, Benin, Burundi, Guinea Bissau, Equatorial Guinea, Mauritius, Eswatini, Djibouti, Reunion, Comoros, Cape Verde, Western Sahara, Tanzania, Togo, Sierra Leone, Libya, Congo, Central African Republic, Liberia, Mozambique, Madagascar, South Africa, Ivory Coast, Cameroon, Niger, Burkina Faso, Mayotte, Sao Tome, Seychelles, South Sudan, Mali, Malawi, Nigeria, Ethiopia, Egypt, DR Congo, Mauritania, Eritrea, Gambia, Botswana, Namibia

#### **Technical Specification :**

DC permanent magnet motor; general notation and mathematical moderlization; Electriocal and mechanical chrastarestics  
Detection of chrastarestic of a mature reaction DC m otor  
Detection cof chrastaretic optoelectronic speed tradecers  
Transducer; potentiometric transducers;  
The study of response of proprtonal?intergrative and derivertive controllers at diffrent input signals



**Laboratory Instrument India**