




Email : sales@laboratoryinstrumentindia.com

Product Name :
Speed and Position Transducer and Control Trainer

Product Code :
LBNY-0006-210001



Description :

Speed and Position Transducer and Control Trainer, Speed and Position Transducer and Control Trainer Equipments, Speed and Position Transducer and Control Trainer Tools, Speed and Position Transducer and Control Trainer Tool Kits, Speed and Position Transducer and Control Trainer Manufacturers, Speed and Position Transducer and Control Trainer Suppliers from India, China, Kenya

Lab Tender Suppliers Africa and Lab Equipment Africa Kenya, Tunisia, Uganda, Sudan, Algeria, Morocco, Angola, Ghana, Gabon, Lesotho, Zambia, Chad, Somalia, Senegal, Zimbabwe, Guinea, Rwanda, Benin, Burundi, Guinea Bissau, Equatorial Guinea, Mauritius, Eswatini, Djibouti, Reunion, Comoros, Cape Verde, Western Sahara, Tanzania, Togo, Sierra Leone, Libya, Congo, Central African Republic, Liberia, Mozambique, Madagascar, South Africa, Ivory Coast, Cameroon, Niger, Burkina Faso, Mayotte, Sao Tome, Seychelles, South Sudan, Mali, Malawi, Nigeria, Ethiopia, Egypt, DR Congo, Mauritania, Eritrea, Gambia, Botswana, Namibia

Technical Specification :

DC permanent magnet motor; general notation and mathematical moderlization; Electriocal and mechanical chrastarestics
Detection of chrastarestic of a mature reaction DC m otor
Detection cof chrastaretic optoelectronic speed tradecers
Transducer; potentiometric transducers;
The study of response of proprtonal?intergrative and derivertive controllers at diffrent input signals



Laboratory Instrument India