




Email : sales@laboratoryinstrumentindia.com

Product Name :
Digital CCTV System

Product Code :
LBNY-0002-670004



Description :

Digital CCTV System

Technical Specification :

Up to 8 Megapixels resolution recording
Support 1-ch HDMI, 1-ch VGA, HMDI at up to 4K(3840x2160) resolution
8-ch network cameras can be connected with 40M/60M/80M incoming bandwidth
Plug & Play with up to 8 independent PoE network interfaces
Support various VCA detection alarm and VCA search
IP Camera Input-8 Supported Cameras
IP Camera up to 8 MP PoE, PoC, Wi-Fi Support: 8xPoE
Internal HDD-2 Alarm input/ Relay Output-4/1 Audio
In/Out-1/1.

Digital CCTV System, Digital CCTV System Equipments, Digital CCTV System Tools, Digital CCTV System Tool Kits, Digital CCTV System Manufacturers, Digital CCTV System Suppliers from India, China, Kenya



Laboratory Instrument India

